

# UNIVERSAL CUTTING MACHINE

Code : B062



- Used to cut construction materials
- It is equipped with a disc holder
- Ideal for wet cutting
- Supplied with precision linear guide bar system with movable specimen trolley, and a disk with dust-proof cover.
- Special custom-made cutting machines (B-062/SP) (with multiple saw discs, custom cutting length ... etc) are also available.

## **SUPPLIED WITH**

- Cutting Blade

## **TECHNICAL SPECIFICATIONS**

- 220 - 240 V / 50 - 60 Hz

## **ORDERING INFORMATION**

<b>Item Name</b>	<b>Item Code</b>
UNIVERSAL CUTTING MACHINE - Ø 350 mm	B062M35XU
UNIVERSAL CUTTING MACHINE - Ø 450 mm	B062M45XU
UNIVERSAL CUTTING MACHINE - Ø 600 mm	B062M60XU
CUTTING BLADE for B-062 - Ø 350 mm	B062P001H
CUTTING BLADE for B-062 - Ø 450 mm	B062P002H
CUTTING BLADE for B-062 - Ø 600 mm	B062P003H
SPARE MOTOR for B062 (3 kW / 1500 rpm)	B062P004U

## **OTHER PHOTOS**

