

UNIVERSAL > PHYSICAL PROPERTIES

# INTERNAL PLASTIC PIPE PRESSURE TESTING MACHINE (Automatic)

Code : U012



*will be added soon ...*

## **TECHNICAL SPECIFICATIONS**

*will be added soon ...*

## **ORDERING INFORMATION**

Item Name	Item Code
INTERNAL PLASTIC PIPE PRESSURE TESTER - 8 bar	U012A008U
PRESSURE CONTROLLING UNIT - 25 bar	U012X025U