

SAMPLE EXTRUDER (Universal, Vertical)

Code : U020 (V)



- Used to extract the sample from its mould. It can be used with proctor test moulds, CBR test moulds, Marshall test moulds and other moulds having the diameter 100 mm to 6".
- The machine is actuated using a manual hydraulic mechanism.

TECHNICAL SPECIFICATIONS

- Manual-Hydraulic Type
- Used for \varnothing 100 mm (\varnothing 4") and \varnothing 150 mm (\varnothing 6") diameter samples
- Can be used with the moulds of Proctor, CBR, Marshall Tests

AVAILABLE MODELS

Code	Item
U 020 H100 H	SAMPLE EXTRUDER, UNIVERSAL - 100 mm
U 020 H300 H	UNIVERSAL SAMPLE EXTRUDER - 300 mm
U 020 H450 H	UNIVERSAL SAMPLE EXTRUDER - 400 mm
U 020 H600 H	UNIVERSAL SAMPLE EXTRUDER - 600 mm
U 020 M300 U	MOTORIZED SAMPLE EXTRUDER - 300 mm

ORDERING INFORMATION

Item Name	Item Code
UNIVERSAL SAMPLE EXTRUDER - 100 mm	U020H100H
UNIVERSAL SAMPLE EXTRUDER - 300 mm	U020H300H
UNIVERSAL SAMPLE EXTRUDER - 450 mm	U020H450H
UNIVERSAL SAMPLE EXTRUDER - 600 mm	U020H600H
MOTORIZED SAMPLE EXTRUDER (Desktop Version) - 100 mm	U020M100U
MOTORIZED SAMPLE EXTRUDER - 300 mm	U020M300U
MOTORIZED SAMPLE EXTRUDER - 650 mm	U020M650U
MOTORIZED SAMPLE EXTRUDER (Horizontal) - 900 mm	U020M900U
ADAPTOR for [U 020] - Ø 38 mm	U020P001H
ADAPTOR for [U 020] - Ø 50 mm	U020P002H
ADAPTOR for [U 020] - Ø 70 mm	U020P003H
ADAPTOR for [U 020] - Ø 106 mm	U020P004H
JACK (30 kN) for U020	U020P005H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 50.80 / 48.50 / 43.50 mm)	U020PH01H

ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 63.50 / 61.00 / 56.00 mm)	U020PH02H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 75.00 / 72.00 / 67.00 mm)	U020PH03H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 83.00 / 79.50 / 74.50 mm)	U020PH04H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 88.90 / 85.00 / 80.00 mm)	U020PH05H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 100.00 / 96.00 / 91.00 mm)	U020PH06H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 101.60 / 97.50 / 92.50 mm)	U020PH07H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 114.30 / 109.00 / 104.00 mm)	U020PH08H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 125.00 / 119.00 / 114.00 mm)	U020PH09H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 106.00 / 100.00 / 95.00 mm)	U020PH10H
ADAPTORS SET for HORIZONTAL EXTRUDER (Ø 156.00 / 150.00 / 145.00 mm)	U020PH11H

OTHER PHOTOS



