

CEMENT TUBE SAMPLER (Packaged)

Code : C028



- The Cement Tube Sampler is specifically designed for obtaining representative samples of hydraulic cement from individual packages, ensuring accurate quality control assessments.
- Constructed from durable brass, the sampler maintains structural integrity and resists corrosion, providing consistent performance across numerous sampling operations.
- Its slender, cylindrical design allows for easy penetration into cement bags, facilitating the extraction of samples without significant disturbance to the remaining material.
- The sampler's dimensions are optimized to collect sufficient quantities of cement for comprehensive testing, aligning with international standards for material sampling.
- Equipped with a comfortable hardwood handle, the sampler offers ease of use during manual operations, reducing operator fatigue and enhancing efficiency.

STANDARDS

ASTM C183 • EN 196-7 • AASHTO T127

TECHNICAL SPECIFICATIONS

- Material: Brass
- Outer Diameter: 32 mm
- Length: 600 mm
- Handle: Hardwood

ORDERING INFORMATION

| Item Name | Item Code |
|--------------------------------|-----------|
| CEMENT TUBE SAMPLER (Packaged) | C028X00XH |