


# LIQUID LIMIT DEVICE (Semi-Automatic, Cone Penetrometer Method)

Code : T090S



*will be added soon ...*

## **TECHNICAL SPECIFICATIONS**

*will be added soon ...*

## **ORDERING INFORMATION**

| Item Name                                | Item Code |
|--|-----------|
| DIGITAL CONE PENETROMETER                | T090H00XB |
| SPARE DIGITAL INDICATOR for PENETROMETER | T090P001H |
| PENETRATION CONE - 30° x 80 g - BS       | T090P002H |

---

|  |           |
|--|-----------|
| PENETRATION CONE - 60° x 60 g - EN       | T090P003H |
| PENETRATION CONE - 30° x 80 g - EN       | T090P004H |
| PENETRATION CONE - 30° x 100 g - EN      | T090P005H |
| PENETRATION CONE - 30° x 400 g - EN      | T090P006H |
| CONE TEST GAUGE - 30°                    | T090P007H |
| CONE TEST GAUGE - 60°                    | T090P008H |
| SLIDING ROD for T090                     | T090P009H |
| SEMI-AUTOMATIC DIGITAL CONE PENETROMETER | T090S00XU |