

WEIGHT SET for CONSOLIDATION and DIRECT SHEAR

Code : T215



- Weight set to be used with consolidation (oedometer) testers and direct shear testers.

TECHNICAL SPECIFICATIONS

- Made of galvanized steel

ORDERING INFORMATION

Item Name	Item Code
SLOTTED WEIGHT for OEDOMETER / DIRECT SHEAR - 1 kg	T215P001H
SLOTTED WEIGHT for OEDOMETER / DIRECT SHEAR - 5 kg	T215P005H

SLOTTED WEIGHT for OEDOMETER / DIRECT SHEAR - 10 kg T215P010H

WEIGHT SET - 50 kg T215X050H

WEIGHT SET - 64 kg T215X064H

WEIGHT SET - 80 kg T215X080H