

# FALLING HEAD PERMEABILITY SET

Code : T250



- Used to determine the permeability of fine-grained soils

## **STANDARDS**

EN 17892-11

## **SUPPLIED WITH**

- Permeameter stand with 3 tubes
- Plated-Steel Falling-head permeability cell - Ø 100 mm (1 ea)
- Water Reservoir Tank

**ORDERING INFORMATION**

Item Name	Item Code
MANOMETER STAND with 4 TUBES (Falling Head)	T250P001H
RESERVOIR TANK (Falling Head)	T250P002H
PERMEABILITY CELL (Falling Head) - Ø 4"	T250P003H
PERMEABILITY CELL (Falling Head) - Ø 6"	T250P004H
MANOMETER TUBE - Ø 6 mm	T250P005H
MANOMETER TUBE - Ø 8 mm	T250P006H
MANOMETER TUBE - Ø 14 mm	T250P007H
MANOMETER TUBE - Ø 22 mm	T250P008H
FALLING HEAD PERMEABILITY SET	T250X00XH