

ASPHALT > MARSHALL STABILITY

MARSHALL COMPACTOR (Automatic, Rotating Mould)

Code : A018



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
MARSHALL COMPACTOR, AUTOMATIC, ROTATING MOULD, 4" - ASTM [50 Hz]	A018M04AT
MARSHALL COMPACTOR, AUTOMATIC, ROTATING MOULD, 4" and 6" - ASTM [50 Hz]	A018M06AT