

## DURIEZ TEST SET

Code : A003



*will be added soon ...*

### **STANDARDS**

*will be added soon ...*

### **TECHNICAL SPECIFICATIONS**

*will be added soon ...*

### **ORDERING INFORMATION**

Item Name	Item Code
-----------	-----------

DURIEZ TESTING MOULD - Ø 80 mm	A003P001H
PENETRATION PISTON - Ø 80 mm	A003P002H
UPPER and LOWER PISTON - Ø 80 mm	A003P003H
TEMP. SUPPORT - Ø 80 mm	A003P004H
DEMOULDING CYLINDRICAL CONTAINER - Ø 80 mm	A003P005H
DURIEZ TESTING MOULD - Ø 120 mm	A003P006H
PENETRATION PISTON - Ø 120 mm	A003P007H
UPPER and LOWER PISTON - Ø 120 mm	A003P008H
TEMP. SUPPORT - Ø 120 mm	A003P009H
DEMOULDING CYLINDRICAL CONTAINER - Ø 120 mm	A003P010H

**OTHER PHOTOS**

