

CENTRIFUGE EXTRACTOR

Code : A031



- Used to determine the bitumen percentage in asphalt mixtures.
- The extractor is available in two capacities (1500 g and 3000 g), and is equipped with a speed control knob.
- *Filter Papers should be ordered separately*

STANDARDS

ASTM D2172 • AASHTO T-164

TECHNICAL SPECIFICATIONS

- Capacity:
 - 1500 g
 - 3000 g
- Speed control knob
- Electronic Brake for fast stopping
- 220 - 240 V / 50 - 60 Hz

ORDERING INFORMATION

Item Name	Item Code
CENTRIFUGE EXTRACTOR - 1500 g	A031M15XU
CENTRIFUGE EXTRACTOR - 3000 g	A031M30XU
FILTER PAPER for A031M15XU	A031P001H
SPARE BOWL for CENTRIFUGE EXTRACTOR - 1500 g	A031P002H
FILTER PAPER for A031M30XU	A031P003H
SPARE BOWL for CENTRIFUGE EXTRACTOR - 3000 g	A031P004H
SHAFT and UPPER BEARING for CENTRIFUGE EXTRACTOR - 3000 g	A031P005H