

# REFLUX EXTRACTOR

Code : A035



- Used to determine the percentage of bitumen in the asphalt pavement mixtures.
- Equipped with a glass cylinder, two wire mesh cones and brass condenser.
- Filter Papers should be ordered separately

## **STANDARDS**

ASTM D2172 • AASHTO T-164

## **TECHNICAL SPECIFICATIONS**

Code	A-035/1	A-035/4
Capacity	1000 g	4000 g
Glass Jar Cylinder	A-035/1/G	A-035/4/G

Wire Mesh Cone (Set of 2)

Thermometer

Hotplate

GL-11/S1

GL-11/S1

Brass Condenser

### **ORDERING INFORMATION**

<b>Item Name</b>	<b>Item Code</b>
CONDENSER for A-035/1	A035P001H
REFLUX BITUMEN EXTRACTOR GLASS - 1000 g	A035P002H
MESH CONES SET for REFLUX BITUMEN EXTRACTOR - 1000 g	A035P003H
CONDENSER for A-035/4	A035P004H
REFLUX BITUMEN EXTRACTOR GLASS - 4000 g	A035P005H
MESH CONES SET for REFLUX BITUMEN EXTRACTOR - 4000 g	A035P006H
FILTER PAPER - Ø 40 cm	A035P007H
REFLUX BITUMEN EXTRACTOR - 1000 g	A035X10XU
REFLUX BITUMEN EXTRACTOR - 4000 g	A035X40XU