

## BITUMEN PENETROMETER (Automatic)

Code : A090A



*will be added soon ...*

### **TECHNICAL SPECIFICATIONS**

*will be added soon ...*

### **ORDERING INFORMATION**

Item Name	Item Code
BITUMEN PENETROMETER - Automatic	A090A00XU
BITUMEN PENETROMETER - Digital	A090H00XB
NEEDLE for BITUMEN PENETROMETER	A090P001H
TRANSFER DISH for PENETROMETER	A090P002H
SLOTTED WEIGHT for PENETROMETER - 50 g	A090P003H

---

SLOTTED WEIGHT for PENETROMETER - 100 g

A090P004H

BITUMEN PENETROMETER - Semi-Automatic

A090S00XU