

FORCE DUCTILITY TESTING MACHINE (Automatic)

Code : A440A



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

| Item Name | Item Code |
|--|-----------|
| FORCE DUCTILITY TESTER | A440A00XU |
| DUCTILITY TESTER with COOLING UNIT [60 Hz] | A440C00XK |
| DUCTILITY TESTER with COOLING UNIT [50 Hz] | A440C00XT |

| | |
|---|-----------|
| DUCTILITY TESTER [60 Hz] | A440M00XK |
| DUCTILITY TESTER [50 Hz] | A440M00XT |
| BASE for DUCTILITY MOULD | A440P001H |
| DUCTILITY MOULD - ASTM | A440P002H |
| ELASTIC RECOVERY DUCTILITY MOULD - EN | A440P003H |
| ELASTIC RECOVERY DUCTILITY MOULD - ASTM | A440P004H |
| FORCE DUCTILITY MOULD - EN | A440P005H |

OTHER PHOTOS

