

# RING AND BALL APPARATUS (Automatic)

Code : A472A



*will be added soon ...*

## **TECHNICAL SPECIFICATIONS**

*will be added soon ...*

## **ORDERING INFORMATION**

Item Name	Item Code
AUTOMATIC RING AND BALL APPARATUS	A472A00XU
RING and BALL SOFTENING POINT APPARATUS - ASTM	A472H00AU
RING and BALL SOFTENING POINT APPARATUS - EN	A472H00EU
STEEL BALL for SOFTENING POINT APPARATUS	A472P001H
BALL CENTERING GUIDE for SOFTENING POINT APPARATUS	A472P002H

---

BRASS RING for SOFTENING POINT APPARATUS

A472P003H