

PELLET PRESS TESTING MACHINE

Code : M002



- Used for forming high-quality pellets for XRF analysis
- Fully automatic with LCD indicator
- Controlling the pressure force, it also determines the time of build-up of force during pressing
- Even difficult materials are pressed perfectly performing pellets with smooth surface
- TFT indicator:
 - Calibration function (password-protected)
 - Maximum load and time are set
 - Overload protection
- Application examples:
 - Cement, gypsum, ores, raw materials, slag ... etc.

TECHNICAL SPECIFICATIONS

- Fully automatic with LCD indicator
- 220 - 240 V / 50 - 60 Hz

ORDERING INFORMATION

Item Name	Item Code
MOULD for M002 - Ø 40 mm	M002P001H
MOULD for M002 - Ø 40 mm	M002P001U
BASE PLATE - Ø 20 cm	M002P002H
MOULD for M002 - Ø 50 mm	M002P003H
MOULD for M002 - Ø 25.23 x 2.70 mm	M002P004H
MOULD for M002 - Ø 14.42 x 3.00 mm	M002P005H
MOULD for M002 - Ø 12.62 x 1.80 mm	M002P006H
PELLET PRESS - 300 kN	M002X00XU
PELLET PRESS - 700 kN	M002X70XU