

# PELLET STRENGTH TESTING MACHINE

Code : M003



- Fully automatic, LCD indicator
- Used to test the strength of pellets
- Can test 60 pcs in a row
- The height is measured before the test
- The average of 60 tests is determined and the test results detects if failure described in the standards
- The graph of each pellet testing is available if required

## **STANDARDS**

- Fully automatic, LCD indicator
- Can test 60 pcs in a row
- 220 - 240 V / 50 - 60 Hz

**ORDERING INFORMATION**

Item Name	Item Code
PELLET STRENGTH TESTER	M003X00XU