

# FLOTATION MACHINE

Code : M050



- Stirrer can be infinitely variable between 50 – 2600 rpm, supplied with digital indicator for speed
- Mechanism to lower/lift the stirring head
- Flotation cells are made of stainless steel with capacities 1, 2, 3 and 4 liters
- The base and support column are made of cast aluminum

## **STANDARDS**

- Adjustable stirring Frequency (50 - 2600 rpm)
- Two flotation stirrers are supplied with the machine:
  - Stainless steel stirrer
  - Polymer stirrer
- Capable of working with 250 – 2000 g of sample
- 220 - 240 V / 50 - 60 Hz

**ORDERING INFORMATION**

<b>Item Name</b>	<b>Item Code</b>
FLOTATION MACHINE	M050X00XU