

## **PRESSURE FILTER**

Code : M053



- Denver-Type Pressure Filter
- Made of stainless steel with a volume of
- Maximum Pressure: 8 kg/cm<sup>2</sup>
- Supplied with air pressure regulator, gauge, valve and filter pack

### **TECHNICAL SPECIFICATIONS**

- Made of stainless steel
- Maximum Pressure: 8 kg/cm<sup>2</sup>

### **ORDERING INFORMATION**

---

<b>Item Name</b>	<b>Item Code</b>
PRESSURE FILTER	M053X00XH