


## **MAGNETIC SEPARATOR - DAVIS-TUBE TYPE**

Code : M074



- Consisting of a powerful electromagnetic field intensity
- The angle of the table is adjustable
- Tube moves forward and backward while rotating
- Supplied with collector

### **TECHNICAL SPECIFICATIONS**

- Wet magnetic separation method
- Adjustable angle

### **ORDERING INFORMATION**

---

Item Name	Item Code
MAGNETIC SEPARATOR - DAVIS-TUBE TYPE	M074X00XU