

# SHREDDER MACHINE

Code : M407



- Used to shred medium-hard minerals, ores, coal, coke, hard plastic, wood, paper, tires and brittle materials
- Heavy-Duty construction
- Low horsepower, energy efficient
- Auto-reversing mechanism to avoid damage in shredding units

## **TECHNICAL SPECIFICATIONS**

- Shredding Discs: 150 - 200 mm in diameter
- Capacity: Depending on material to be shredded, 50 - 250 kg/hr
- 380 V / 50 Hz

**ORDERING INFORMATION**

<b>Item Name</b>	<b>Item Code</b>
SHREDDER MACHINE	M407X00XV