

## BALL MILL

Code : M438



- Used to dry-grind lime stones, ceramic, glass and slag

### TECHNICAL SPECIFICATIONS

- Drum made of stainless steel (300 x 300 mm)
- Adjustable Rotation Speed (Max. 70 rpm)
- Supplied with:
- Grinding Balls:
  - 43 x Ø38.1 mm
  - 67 x Ø31.75 mm
  - 10 x Ø25.4 mm
  - 71 x Ø19.05 mm
  - 94 x Ø15.87 mm

- Grinding Mills:
  - 10 x 20 mm
  - 10 x 25 mm
  - 10 x 30 mm
- 220-240 V/ 50-60 Hz

**ORDERING INFORMATION**

Item Name	Item Code
ABRASIVE CHARGES for M438	M438P001H
ABRASIVE CHARGES for M438 - Ø 50 mm / Pack of 5 kg	M438P002H
ABRASIVE CHARGES for M438 - Ø 40 mm / Pack of 5 kg	M438P003H
ABRASIVE CHARGES for M438 - Ø 30 mm / Pack of 5 kg	M438P004H
ABRASIVE CHARGES for M438 - Ø 25 mm / Pack of 1 kg	M438P005H
ABRASIVE CHARGES for M438 - Ø 19 mm / Pack of 1 kg	M438P006H
ABRASIVE CHARGES for M438 - Ø 16 mm / Pack of 1 kg	M438P007H
ABRASIVE CHARGES for M438 - Ø 12.5 mm / Pack of 1 kg	M438P008H
BALL MILL	M438X00XU