

# HAMMER CRUSHER

Code : M440



- Used to crush medium-hard materials and specially coal and coke
- Supplied with feed hopper and collector
- Automatic cut-off electricity if the
- Sieves (3 ... 9 mm) upon request (to be specified at the time of order)
- Easy to clean

## **TECHNICAL SPECIFICATIONS**

- Material feed size is < 50 mm
- Final fineness is maximum 3 mm (with 3 mm sieve installed)
- Capacity: 250 – 700 kg/hr
- 220 - 240 V / 50 - 60 Hz

**ORDERING INFORMATION**

<b>Item Name</b>	<b>Item Code</b>
HAMMER CRUSHER	M440X00XV