

# VACUUM PUMP

Code : G702



Used to vacuum specimens, remove air voids and speed up deairing process.

## **ORDERING INFORMATION**

Item Name	Item Code
HOSE for VACUUM PUMP	G702P001H
GAUGE for VACUUM PUMP	G702P002H
REGULATOR for VACUUM PUMP	G702P003H
VACUUM PUMP - 51 lt/min	G702X051U
VACUUM PUMP - 100 lt/min	G702X100U
VACUUM PUMP - 128 lt/min (Dual Stage)	G702X128U
VACUUM PUMP - 170 lt/min (Dual Stage)	G702X170U