

## HAIR BRUSH

Code : G719



- A hair brush in engineering testing laboratories is used for gently cleaning fine mesh sieves without damaging them. It helps remove fine particles stuck in the mesh, ensuring accurate sieve analysis and maintaining the integrity of the sieve.

### **ORDERING INFORMATION**

Item Name	Item Code
HAIR BRUSH - No: 1	G719X001H
HAIR BRUSH - No: 2	G719X002H
HAIR BRUSH - No: 3	G719X003H
HAIR BRUSH - No: 4	G719X004H
BRUSH with HANDLE	G719X005H