

## MOISTURE TIN

Code : G720



- Moisture tins in engineering testing laboratories are used to hold soil or material samples during oven drying for moisture content determination. They are durable, heat-resistant, and help ensure accurate weight measurements before and after drying.

### **ORDERING INFORMATION**

Item Name	Item Code
MOISTURE TIN, with LID - Small - Ø 55 x 35 mm	G720X001H
MOISTURE TIN, with LID - Medium	G720X002H
MOISTURE TIN, with LID - Large	G720X003H
MOISTURE TIN - Ø 200 x 50 mm	G720X011H
MOISTURE TIN - Ø 75 x 30 mm	G720X012H