

# WATER QUALITY TEST KIT

Code : G760



- A water quality test kit in engineering testing laboratories is used to assess various parameters of water, such as pH, turbidity, dissolved oxygen, and chemical contaminants. It helps ensure that water meets environmental and regulatory standards for quality before use in construction or other processes.

## **ORDERING INFORMATION**

Item Name	Item Code
WATER QUALITY TEST KIT	G760X001H
WATER QUALITY TESTING SET	G760X002B
SALT RATIO MEASUREMENT KIT	G760X002H
SALT RATIO MEASUREMENT DEVICE	G760X003B