

## SPIRIT LEVEL

Code : G764



- A spirit level in engineering testing laboratories is used to check the horizontal or vertical alignment of surfaces or equipment. It ensures accurate positioning during sample preparation, equipment setup, and testing procedures, improving the reliability of results.

### **ORDERING INFORMATION**

Item Name	Item Code
SPIRIT LEVEL - 20 cm	G764X020H
SPIRIT LEVEL - 30 cm	G764X030H
SPIRIT LEVEL - 60 cm	G764X060H
SPIRIT LEVEL - 80 cm	G764X080H
SPIRIT LEVEL - 120 cm	G764X120H

---

SPIRIT LEVEL - 180 cm

G764X180H