

## SAMPLE TRAY

Code : G768



- A sample tray in engineering testing laboratories is used to hold, transport, and organize material samples during testing or analysis. It helps prevent contamination, ensures proper handling, and allows for easy access and examination of samples.

### **ORDERING INFORMATION**

Item Name	Item Code
SAMPLE TRAY - X Small Size, 200 x 300 x 40 mm	G768X001H
SAMPLE TRAY - Small Size, 230 x 330 x 40 mm	G768X002H
SAMPLE TRAY - Medium Size, 260 x 380 x 40 mm	G768X003H
SAMPLE TRAY - Large Size, 280 x 430 x 40 mm	G768X004H
SAMPLE TRAY - X Large Size, 300 x 480 x 40 mm	G768X005H

---

TRAY -

G768X00XH