

# PLASTIC CONTAINER

Code : W152



- A plastic container in engineering testing laboratories is used for storing, transporting, and organizing samples, liquids, chemicals, or materials. It is lightweight, durable, and resistant to breakage, providing a safe and practical solution for holding substances during experiments or when not in use.

## **ORDERING INFORMATION**

Item Name	Item Code
PLASTIC CONTAINER - 5 lt	W152P005H
PLASTIC BUCKET - 20 lt	W152P020H
PLASTIC BUCKET - 30 lt	W152P030H