

TRICHLOROETHYLENE

Code : G801



- Trichloroethylene (TCE) in engineering testing laboratories is commonly used as a solvent for cleaning equipment, dissolving fats, oils, and greases, and as a degreasing agent in various industrial applications. It can also be used in chemical analysis and testing, although it requires careful handling due to its toxic and carcinogenic properties.

ORDERING INFORMATION

Item Name	Item Code
TRICHLOROETHYLENE	G801X00XH
TRICHLOROETHYLENE - 1 Barrel	G801X01XH