

QUANTAB CHLORIDE TITRATOR STRIPS (Type 1175)

Code : G815



- The Quantab Chloride Titrator 1175 type in engineering testing laboratories is used for determining the chloride content in water, soil, or concrete samples. It is a quick and efficient method for assessing the chloride levels, which is crucial for evaluating corrosion potential, especially in materials like reinforced concrete.

ORDERING INFORMATION

Item Name	Item Code
QUANTAB CHLORIDE TITRATOR 1175 TYPE	G815X00XH