

QUANTAB CHLORIDE TITRATOR STRIPS (Type 1176)

Code : G816



- The Quantab Chloride Titrator 1176 type in engineering testing laboratories is used to measure chloride concentrations in water, concrete, and soil samples. It is designed for accurate, quick, and reliable titration, helping assess the risk of corrosion in materials, particularly in concrete structures exposed to chloride contamination.

ORDERING INFORMATION

Item Name	Item Code
QUANTAB CHLORIDE TITRATOR 1176 TYPE	G816X00XH