

# VOLUMETRIC FLASK

Code : W003



- A volumetric flask in engineering testing laboratories is used for accurately preparing and diluting solutions to a specific volume. Its precision makes it ideal for measuring exact concentrations in chemical, environmental, or material testing, ensuring consistency and accuracy in experiments.

## **ORDERING INFORMATION**

Item Name	Item Code
VOLUMETRIC FLASK - 5 ml	W003G005H
VOLUMETRIC FLASK - 10 ml	W003G010H
VOLUMETRIC FLASK - 1000 ml	W003G01KH
VOLUMETRIC FLASK - 25 ml	W003G025H

VOLUMETRIC FLASK - 2000 ml	W003G02KH
VOLUMETRIC FLASK - 3000 ml	W003G03KH
VOLUMETRIC FLASK - 4000 ml	W003G04KH
VOLUMETRIC FLASK - 50 ml	W003G050H
VOLUMETRIC FLASK - 5000 ml	W003G05KH
VOLUMETRIC FLASK - 6000 ml	W003G06KH
VOLUMETRIC FLASK - 100 ml	W003G100H
VOLUMETRIC FLASK - 150 ml	W003G150H
VOLUMETRIC FLASK - 250 ml	W003G250H
VOLUMETRIC FLASK - 500 ml	W003G500H

**OTHER PHOTOS**

