

# LABORATORY BOTTLE

Code : W006



- A laboratory bottle in engineering testing laboratories is used for storing, transporting, and handling liquids, chemicals, or samples. These bottles are typically made of glass or durable plastic and are designed to safely contain substances while providing ease of access and measurement during testing procedures.

## **ORDERING INFORMATION**

Item Name	Item Code
LABORATORY BOTTLE - 1000 ml	W006G01KH
LABORATORY BOTTLE - 25 ml	W006G025H
LABORATORY BOTTLE - 2000 ml	W006G02KH

LABORATORY BOTTLE - 50 ml	W006G050H
LABORATORY BOTTLE - 100 ml	W006G100H
LABORATORY BOTTLE - 250 ml	W006G250H
LABORATORY BOTTLE - 500 ml	W006G500H

**OTHER PHOTOS**

