

PETRI DISH

Code : W007



- A Petri dish in engineering testing laboratories is used to culture microorganisms or hold samples for observation. It is commonly used in microbiological studies, material testing, and surface analysis, providing a flat, sterile surface for samples or experiments.

ORDERING INFORMATION

Item Name	Item Code
PETRI DISH - 60 mm	W007G060H
PETRI DISH - 80 mm	W007G080H
PETRI DISH - 100 mm	W007G100H
PETRI DISH - 120 mm	W007G120H
PETRI DISH - 150 mm	W007G150H