

VACUUM DESSICATOR

Code : W009



- A vacuum desiccator in engineering testing laboratories is used to store samples in a low-pressure, dry environment. By applying a vacuum, it removes moisture from the air and from the sample itself, ensuring that samples remain dry and uncontaminated during storage or preparation for sensitive tests.

ORDERING INFORMATION

Item Name	Item Code
VACUUM DESSICATOR - 200 mm	W009G296H
VACUUM DESSICATOR - 250 mm	W009G344H
VACUUM DESSICATOR - 300 mm	W009G420H