


GLASS MEASURING CYLINDER

Code : W010



- A glass measuring cylinder in engineering testing laboratories is used for accurately measuring the volume of liquids or small amounts of substances. Its clear markings and wide cylindrical shape ensure precise volume measurements during experiments, chemical analysis, or sample preparation.

ORDERING INFORMATION

Item Name	Item Code
GLASS MEASURING CYLINDER - 10 ml	W010G010H
GLASS MEASURING CYLINDER - 1000 ml	W010G01KH
GLASS MEASURING CYLINDER - 25 ml	W010G025H
GLASS MEASURING CYLINDER - 2000 ml	W010G02KH

GLASS MEASURING CYLINDER - 50 ml	W010G050H
----------------------------------	-----------

GLASS MEASURING CYLINDER - 100 ml	W010G100H
-----------------------------------	-----------

GLASS MEASURING CYLINDER - 250 ml	W010G250H
-----------------------------------	-----------

GLASS MEASURING CYLINDER - 500 ml	W010G500H
-----------------------------------	-----------

OTHER PHOTOS

