

## PYCNOMETER (Gay Lussac)

Code : W014



- A pycnometer (Gay Lussac type) in engineering testing laboratories is used to determine the specific gravity or density of liquids or powders. It consists of a precisely calibrated glass container, often with a stopper, ensuring accurate measurements by accounting for temperature variations during testing. This is essential for density and quality control in chemical, material, and environmental testing.

### **ORDERING INFORMATION**

Item Name	Item Code
PYCNOMETER (Gay Lussac) - 5 ml	W014G005H
PYCNOMETER (Gay Lussac) - 10 ml	W014G010H
PYCNOMETER (Gay Lussac) - 1000 ml	W014G01KH

PYCNOMETER (Gay Lussac) - 25 ml	W014G025H
PYCNOMETER (Gay Lussac) - 2000 ml	W014G02KH
PYCNOMETER (Gay Lussac) - 50 ml	W014G050H
PYCNOMETER (Gay Lussac) - 100 ml	W014G100H
PYCNOMETER (Gay Lussac) - 250 ml	W014G250H
PYCNOMETER (Gay Lussac) - 500 ml	W014G500H

**OTHER PHOTOS**

