


PLASTIC MEASURING CYLINDER

Code : W032



- A plastic measuring cylinder in engineering testing laboratories is used for accurately measuring the volume of liquids. It is lightweight and durable, making it suitable for handling a variety of chemicals or materials, especially when glass containers might be too fragile or heavy. It provides precise measurements for experiments and sample preparation.

ORDERING INFORMATION

Item Name	Item Code
PLASTIC MEASURING CYLINDER - 10 ml	W032P010H
PLASTIC MEASURING CYLINDER - 1000 ml	W032P01KH
PLASTIC MEASURING CYLINDER - 25 ml	W032P025H
PLASTIC MEASURING CYLINDER - 2000 ml	W032P02KH

PLASTIC MEASURING CYLINDER - 50 ml	W032P050H
------------------------------------	-----------

PLASTIC MEASURING CYLINDER - 100 ml	W032P100H
-------------------------------------	-----------

PLASTIC MEASURING CYLINDER - 250 ml	W032P250H
-------------------------------------	-----------

PLASTIC MEASURING CYLINDER - 500 ml	W032P500H
-------------------------------------	-----------