

## TEST TUBE

Code : W038



*will be added soon ...*

### **TECHNICAL SPECIFICATIONS**

*will be added soon ...*

### **ORDERING INFORMATION**

Item Name	Item Code
TEST TUBE (Height: 75 mm)	W038G075H
TEST TUBE (Height: 100 mm)	W038G100H
TEST TUBE (Height: 160 mm)	W038G160H
TEST TUBE (Height: 180 mm)	W038G180H
TEST TUBE (Height: 200 mm)	W038G200H