

MOBILE LABORATORY

Code : P001



- We can offer a complete line of Mobile Laboratory with special design to meet all possible environments such as cold or tropical climate.
- The main tree types are:
 - Container
 - Trailer Mounted
 - Van Mounted

TECHNICAL SPECIFICATIONS

- Office / Bedroom section
- Lighting & Electric sockets
- Electric generator
- Air Conditioning / Aspirator Fan

- Water Tank, WC & Sink

ORDERING INFORMATION

Item Name	Item Code
TRAILER for MOBILE LABORATORY - 20 DC	P001P001X
MOBILE LABORATORY - 20 DC	P001X00XX

OTHER PHOTOS

