

ASPHALT > MARSHALL STABILITY

MARSHALL COMPACTOR (Automatic, Triple)

Code : A017



- Used to determine the resistance of asphalt sample to plastic flow. The test is performed in a fully automatic mechanism and programmed to stop automatically at the required number of blows.
- Supplied with the standard height and rammer weight according to ASTM/EN standards.
- Equipped with a blow counter, wooden pedestal and quick action clamping.

STANDARDS

ASTM D6926 • EN 12697-10 • AASHTO t245

TECHNICAL SPECIFICATIONS

- Consists of 3-mould structure that can compact 3 specimens at the same time
- Blow Rate: 55 ± 5 blows per minute

- Rammer Diameter: 100.33 mm $\pm$ 0.25 mm
- Rammer Weight: 4536 g $\pm$ 9 g
- Drop Height: 457.205 mm $\pm$ 0.635 mm
- 220 - 240 V / 50 - 60 Hz

ORDERING INFORMATION

Item Name	Item Code
TRIPLE MARSHALL COMPACTOR, AUTOMATIC, STATIONARY BASE - ASTM [60 Hz]	A017M04AK
TRIPLE MARSHALL COMPACTOR, AUTOMATIC, STATIONARY BASE - ASTM [50 Hz]	A017M04AT