

SAYBOLT VISCOMETER

Code : A430



- Used to determine the viscosity of the petroleum products at a specified temperature.
- The body of the machine is made of stainless steel.

STANDARDS

ASTM D88 • AASHTO T-72

EQUIPPED WITH

- Electrical Heater with Digital Thermostat
- Mixer for Homogeneity

SUPPLIED WITH

- Filter funnel
- Universal Orifice
- Furol Orifice
- Saybolt receiving flask

ORDERING INFORMATION

Item Name	Item Code
UNIVERSAL ORIFICE for SAYBOLT VISCOMETER	A430P002H
FUROL ORIFICE for SAYBOLT VISCOMETER	A430P003H
FILTER FUNNEL w/ WIRE MESH for SAYBOLT VISCOMETER	A430P004H
WITHDRAWAL TUBE for SAYBOLT VISCOMETER	A430P005H
CALIBRATION OIL for SAYBOLT - 500 ml	A430P006H
CALIBRATION OIL for SAYBOLT - 60 ml	A430P007H
SAYBOLT VISCOMETER	A430X01AU
SAYBOLT VISCOMETER - Two-Tube	A430X02AU