

ASPHALT > BITUMEN PROPERTIES

SOLUBILITY TESTING MACHINE

Code : A438



- Used to determine the degree of solubility in trichloroethylene of asphalt materials having little or no mineral matter.

STANDARDS

ASTM D2042 • EN 12592

SET

- Filter Flask - 500 ml with Rubber Stopper
- Erlenmeyer - 500 ml
- Vacuum Pump with 1 m flexible hose
- Gooch Crucible with Funnel

-
- Filter Papers (Ø 25 mm) - 5 x Pack of 100

ORDERING INFORMATION

Item Name	Item Code
SOLUBILITY TESTER	A438X00XU