

ASPHALT > DUCTILITY

FORCE DUCTILITY TESTING MACHINE (Automatic)

Code : A440A



will be added soon ...

TECHNICAL SPECIFICATIONS

will be added soon ...

ORDERING INFORMATION

Item Name	Item Code
FORCE DUCTILITY TESTER	A440A00XU
DUCTILITY TESTER with COOLING UNIT [60 Hz]	A440C00XK
DUCTILITY TESTER with COOLING UNIT [50 Hz]	A440C00XT

DUCTILITY TESTER [60 Hz]	A440M00XK
DUCTILITY TESTER [50 Hz]	A440M00XT
BASE for DUCTILITY MOULD	A440P001H
DUCTILITY MOULD - ASTM	A440P002H
ELASTIC RECOVERY DUCTILITY MOULD - EN	A440P003H
ELASTIC RECOVERY DUCTILITY MOULD - ASTM	A440P004H
FORCE DUCTILITY MOULD - EN	A440P005H

OTHER PHOTOS

