

ASPHALT > BITUMEN PROPERTIES

FLASH and FIRE POINT APPARATUS (Cleveland Open Cup)

Code : A470



- Used to determine the maximum temperature bitumen can be safely heated without the danger of an instantaneous flash in the presence of an open flame.
- Supplied with an electronic heater, a thermoregulator, two thermometers, and a sample cup.

STANDARDS

ASTM D117 • ASTM D92 • AASHTO T-48

SUPPLIED WITH

- Electric heater

- Thermoregulator
- Two Thermometers
- Sample Cup

ORDERING INFORMATION

Item Name	Item Code
AUTOMATIC FLASH and FIRE POINT APPARATUS (Cleveland Open Cup)	A470A00XU
FLASH and FIRE POINT APPARATUS (Cleveland Open Cup)	A470H00XU