

FLEXURAL TESTING MACHINE (C-Frame)

Code : B002



- Used to determine the flexural strength of the test specimen
- Different models are available:
 - C-Frame
 - Load Sensor: Loadcell
 - Rectangular Frame
 - Load Sensor: Pressure Transducer
- Different Capacities are available:
 - 200 kN
 - 300 kN
 - 600 kN
- Can be purchased as stand alone machine or connected to an existing hydraulic unit (must be mentioned at the time of order)

ORDERING INFORMATION

Item Name	Item Code
FLEXURAL TESTING MACHINE (C-Type) - 200 kN	B002A20XU
FLEXURAL TESTING MACHINE (C-Type) - 300 kN	B002A30XU
POWER PACK (for Compression Frame) - 400 Bar	B002P000U
POWER PACK (for Flexure Frame) - 400 Bar	B002P001U
POWER PACK (for Compression Frame) - 400 Bar (Without Controlling Unit)	B002P002U
POWER PACK (Servo-Motor, for Compression Frame) - 400 Bar	B002P003U
HYDRAULIC UNIT (for Compression Frame) - 400 Bar	B002P004U
POWER PACK (for Semi-Automatic Compression Frame) - 400 Bar	B002P005U
POWER PACK with PILOT VALVE - 400 Bar (Without Controlling Unit)	B002P006U
UPGRADE PACK for MODULUS OF ELASTICITY MEASUREMENT	B002P007H
POWER PACK with PILOT VALVE - 400 Bar	B002P008U
FLEXURAL FRAME (C-Type) - 200 kN	B002P20XH
FLEXURAL FRAME (C-Type) - 300 kN	B002P30XH

OTHER PHOTOS

