

CYLINDER MOULD

Code : B010



- Available as:
 - Plastic
 - Steel
- Different dimensions are available

ORDERING INFORMATION

Item Name	Item Code
CYLINDER MOULD - Metal (Ø 10 cm)	B010P001H
CYLINDER MOULD - Metal (Ø 15 cm)	B010P002H
CYLINDER MOULD - Metal (Ø 16 cm, Thickness: 4 mm)	B010P003H

CYLINDER MOULD - Metal (Ø 16 cm, Thickness: 8 mm)	B010P004H
CYLINDER MOULD - Metal (Ø 25 cm)	B010P005H
CYLINDER MOULD - Plastic (Ø 10 cm)	B010P006H
CYLINDER MOULD - Plastic (Ø 15 cm)	B010P007H
CYLINDER MOULD - Plastic (Ø 16 cm)	B010P008H
CYLINDER MOULD - Metal (Ø 6" x 12" / Ø 152 x 305 mm)	B010P009H
PLUG for PLASTIC CYLINDER MOULD	B010P010H
CYLINDER MOULD (SPRING LID TYPE) - 100 x 200 mm	B010P015H
CYLINDER MOULD (BILLY LID TYPE) - 100 x 200 mm	B010P025H

OTHER PHOTOS

