

CONCRETE > FRESH CONCRETE TESTS

KELLY-BALL PENETRATION APPARATUS

Code : B022



- Used to determine the workability of concrete as an alternative to Slump Test, where 1" of kelly ball penetration corresponds to 2" of slump value
- Coated for rust protection
- Weight of the movable parts: 30 lb (13.6 kg)

STANDARDS

ASTM C360

TECHNICAL SPECIFICATIONS

- Metal hemispherical cylinder:
 - Diameter : Ø 6"

- Height : 4 5/8"

SUPPLIED WITH

- Metal hemispherical cylinder and graduated handle

ORDERING INFORMATION

Item Name	Item Code
KELLY BALL APPARATUS	B022X00AH