

CONCRETE > SELF-CONSOLIDATING CONCRETE

## CONCRETE BLEEDING TEST SET

Code : B028



- The Concrete Bleeding Test Set, compliant with EN 480-4 standards, is essential for assessing the bleeding characteristics of freshly mixed concrete.
- Bleeding refers to the migration of water to the surface of concrete, which can impact its quality, durability, and workability.

### **STANDARDS**

EN 480-4

### **TECHNICAL SPECIFICATIONS**

- Cylindrical container with Cover (Ø 250 x 280 mm)
  - Tamping Rod
  - Pipette (50 ml)
  - Glass Beaker (100 ml)

**ORDERING INFORMATION**

| Item Name                  | Item Code |
|----------------------------|-----------|
| CONCRETE BLEEDING TEST SET | B028X00EH |